

CATALOGUE | PASBOXX

PASSBOXX NDS 1600

NADISOLAB

Pass box is a cubicle box which has got interlock doors located on both sides



PT. KARYA NADISO UTAMA

www.nadisolab.com

NDS

I-600

ABOUT PASS BOX NDS I-600


As a kind of auxiliary equipment used in clean room or sterilize room, pass box is mainly used to transfer small object (media sample) between common area and clean room,

Pass Box Nadiso used interlock to open one side and lock other side, and it operate by system touch screen with dimension 2,8".



WHY CHOOSE PASS BOX NDS ?

- Our pass boxes have an ergonomic shape and design
- touch screen system makes operation easy
- there is UV which can minimize the development of bacteria
- the position of the door that is parallel to the body, this makes the room design so spacious





TOP VIEW

SPECIFICATION PASS BOX

ISOMET

Nadiso Pass Box is a cubicle box which has got interlocked doors located on both sides. This protects the ® controlled environment from being polluted while the transfer of a material is taking place inside. Pass boxes are the perfect tool for transporting items in and out of cleanrooms. These boxes are widely used in laboratories, hospitals, pharmaceutical and light-industrial environments. They not only offer full protection but also facilitate transportation of products.

- External Size
- Shower Area
- Material Specification
- Display
- Interlock
- UV Lamp
- HEPA Filter
- Clean Air Speed
- Fluorescent
- Power Supply
- 600(L) x 600(W) x 580(H) mm
- 500(L) x 500(W) x 500(H) mm
- Full Stainless Steel 1,2 & 2 mm Thick
- LCD Touch Screen 2,8 inch
- Magnet Interlock system
- 30 Watt
- No
- No
- No
- 220V, 60Hz

500


SIDE VIEW



FRONT VIEW

Catalogue
Nadisolab

CATALOGUE PASS BOX NDS I-600



For more information, please contact
PT. Karya Nadiso Utama

Head Office & Wokshop

Kawasan Industri Manis

Jl. Manis IV No.5 RT 01/03 Kel. Manis Jaya

Kec Jatiuwung, Tangerang 15136

Telp : +62 (21) 5421 0226 - 28

Fax : +62 (21) 5421 0229

WareHouse

Jl. Diklat Pemda, Dukuh Pinang No.24

RT.001/003 Kel. Bojong Nangka

Kec Kelapa Dua Tangerang 15821